OVITEX® CORE

REINFORCED TISSUE MATRIX

POLYMER REINFORCEMENT
RESORBABLE Poly(glycolic Acid)

Distributed by:

1 Great Valley Parkway, Suite 24
Malvern, PA 19355 USA
T (844) 835-2246
F (844) 455-8707
www.telabio.com

Aroa Biosurgery Ltd.
2 Kingsford Smith Place,
Airport Oaks, Auckland 2022,
New Zealand

S4M EUROPE
59, rue Castellion
01100 Oyonnax - France

AROA™ and AROA ECM™ are trademarks of Aroa Biosurgery Limited.
TELA Bio and OviTex® are registered trademarks
of TELA Bio, Inc. © TELA Bio, Inc.

OVITEX® CORE

REINFORCED TISSUE MATRIX

POLYMER REINFORCEMENT
RESORBABLE Poly(glycolic Acid)

SCIENCE. VALUE. INNOVATION.

OviTex® Reinforced Tissue Matrix

Resorbable Polymer

EN

INSTRUCTIONS FOR USE

DESCRIPTION OF THE PRODUCT

OviTex® Core Reinforced Tissue Matrix with Resorbable Polymer (OviTex Core) is a sterile surgical mesh composed of ovine (sheep) derived extracellular matrix (ECM) and polyglycolic acid (PGA). OviTex Core will incorporate into the recipient tissue with associated cellular and microvascular ingrowth.

OviTex Core is provided in various sizes to suit surgeon preference and the complexity of the soft tissue repair. The rectangular device may be trimmed to a desired shape to further accommodate an individual patient's requirements.

INDICATIONS FOR USE

OviTex Core is intended for use as a surgical mesh to reinforce and/or repair soft tissue where weakness exists. Indications for use include the repair of hernias and/or abdominal wall defects that require the use of reinforcing or bridging material to obtain the desired surgical outcome.

CONTRAINDICATIONS

Do not use OviTex Core in patients with a known sensitivity to materials of ovine (sheep) origin. Use of OviTex Core in this patient population may result in an allergic or immunological reaction.

RISK STATEMENTS

WARNINGS

- Device is supplied sterile. Inspect the packaging to ensure it is intact and undamaged prior to use. Do not use OviTex Core, and discard prior to use if the packaging is damaged. Do not re-sterilize.
- Single use only. Reuse, resterilization, reprocessing and/or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the device, which may result in device failure and/or patient injury. Open and unused material should be discarded.

PRECAUTIONS

- Do not use the product past its expiration date. The expiration date is displayed on the product labeling as the year (4 digits), month (2 digits), and day (2 digits) next to an hourglass symbol.
- Place the device in maximum contact with healthy, well-vascularized tissue to promote cell ingrowth and tissue remodeling.

ADVERSE EVENTS

The following adverse events have been reported for surgical repair of hernias (with or without a surgical mesh): pain, infection, hernia recurrence, adhesion, bowel obstruction, bleeding, fistula, seroma, perforation, mesh migration, and mesh contraction.

MRI SAFETY INFORMATION

OviTex Core is MR Safe.

HOW SUPPLIED

OviTex Core is packaged in a Tyvek/film double pouch configuration.

LATEX INFORMATION

OviTex Core is not made with natural rubber latex.

STORAGE

OviTex Core should be stored at room temperature in a clean and dry area.

CAUTION

OviTex Core is available by prescription only.

INSTRUCTIONS

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

1. Inspect the packaging to ensure it is intact and undamaged.
2. Using aseptic technique, remove the inner pouch from its outer pouch and place the inner pouch in the sterile field.
3. Open the inner pouch carefully and aseptically remove the device using sterile forceps.
4. Place the device into a sterile dish in the sterile field.
5. Rehydrate the device in a sufficient volume of sterile saline or sterile Lactated Ringer's solution for a minimum of 5 minutes.
6. Prepare the site using standard surgical techniques.
7. Using aseptic technique, trim the device to fit the site, if necessary, providing an allowance for overlap. Position the device to achieve maximum contact between the device and surrounding tissue. To facilitate cell migration and tissue ingrowth, an overlap of 3-5 cm with healthy well-vascularized tissue is suggested.

Note: If the device is cut too small for the defect, excess tension may be placed on the suture line. This can result in recurrence of the original tissue defect or development of a defect in the adjacent tissue.

8. Using aseptic technique, transfer the device to the surgical site and suture, staple, or tack into place, avoiding excess tension.
9. Complete the surgical procedure.
10. Discard any unused portions according to institutional guidelines for medical waste.

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SYMBOLS GLOSSARY

Symbols contained in ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements are indicated below.		
	5.1.1	Manufacturer
	5.1.2	European Commission Authorized Representative
	5.1.3	Date of manufacture
	5.1.4	Use by/Expiration date
	5.1.5	Lot number
	5.1.6	Catalog number
	5.1.8	Importer
	5.2.3	Sterilized using ethylene oxide
	5.2.6	Do not resterilize
	5.2.8	Do not use if packaging is damaged
	5.2.12	Double sterile barrier system
	5.3.4	Keep dry
	5.4.2	Do not reuse
	5.4.3	Consult instructions for use
	5.4.8	Contains biological material of animal origin
	5.7.7	Medical Device
	-	Number of units
Rx Only	-	Prescription use only. Federal (USA) law restricts this device to sale by or on the order of a physician.

EN
Refer to Eudamed for a summary of safety and clinical performance.
Report any serious incident related to OviTex Core to TELA Bio via S4M Europe, as well as the applicable Member State competent authority.

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